


ADIAVET™ CEMO REAL TIME

QUICK DETECTION OF CONTAGIOUS EQUINE METRITIS BY REAL-TIME PCR

ADIA GENE

Contagious equine metritis is a sexually transmitted disease specific to equids caused by the bacterium *Taylorella equigenitalis* that is notifiable to the OIE. This disease can have reproductive consequences (sterility and sometimes abortions). Infected stallions are asymptomatic. In mares the infection is characterized by abundant vaginal exudate. Other bacteria, such as *Taylorella asinigenitalis*, *Klebsiella pneumoniae*, and *Pseudomonas aeruginosa*, are also found in metritis and may require treatment.

The combination of a **quick protocol** for treating the samples (swabs, semen, bacterial cultures) and two triplex real-time PCR kits, **ADIAVET™ CEMO TAYLORELLA REAL TIME** and **ADIAVET™ CEMO KLEB/PSEUDO REAL TIME**, allows rapid, specific, sensitive identification of these four pathogens.



ADVANTAGES

- Sample preparation method that is quick, easy, and inexpensive
- Ready-to-use PCR reagents
- Kits stable for 12 months at -20 °C
- DNA stable for 24 hours at 4 °C
- Swabs stable for 7 days at room temperature
- Possibility of using swabs without charcoal friends
- Possibility of extraction on silica columns
- Validated according to standard NF U47-600

TRIPLEX PCR KIT REFERENCES

	ADIAVET™ CEMO TAYLORELLA REAL TIME (réf. ADI501-100)	ADIAVET™ CEMO KLEB/PSEUDO REAL TIME (Réf. ADI511-100)
FAM	<i>Taylorella equigenitalis</i>	<i>Klebsiella pneumoniae</i>
CY5	<i>Taylorella asinigenitalis</i>	<i>Pseudomonas aeruginosa</i>
HEX	IPC amplification	IPC amplification


- Positive controls included in kits
- DL_{PCR} and sentinel controls for help in adopting the method available on request

PERFORMANCE

PCR

- 100% inclusiveness (more than 20 strains/targets tested)
- 100% exclusiveness (more than 100 strains tested)
- $DL_{PCR} \leq 20$ GE/PCR GE: genome equivalent

DL OF THE METHOD



DIAGNOSTIC SENSITIVITY AND SPECIFICITY

500 genital swabs tested by culturing (AFNOR U47-108) and by PCR with the direct lysis method at LABÉO Frank Duncombe

	Se	Sp
<i>K. pneumoniae</i>	100%	93%
<i>P. aeruginosa</i>	100%	99%
<i>T. equigenitalis</i>	Not determined*	100%
<i>T. asinigenitalis</i>	Not determined*	100%


* No *T. equigenitalis*- and *T. asinigenitalis*-positive samples in the entire study.

VALIDATIONS

> Validation of *Taylorella* real-time PCR on genital swabs

Negative genital swabs doped with *T. equigenitalis* or *T. asinigenitalis* strains

Blind assessments at AdiaGène and LABÉO Frank Duncombe



- 100% of positives detected
- No false positives
- Direct lysis and purification methods give equivalent results
- Good repeatability and reproducibility (CV < 3%)

> Validation of two real-time PCR kits on qualified samples

PCR testing of 9 samples from an inter-laboratory culture trial (Horserace Betting Levy Board)

Culture	ADIAVET™ CEMO TAYLORELLA REAL TIME		ADIAVET™ CEMO KLEB/PSEUDO REAL TIME	
	<i>T. equigenitalis</i>	<i>T. asinigenitalis</i>	<i>K. pneumoniae</i>	<i>P. aeruginosa</i>
#1 <i>P. aeruginosa</i>	-	-	-	+ (Ct = 26,5)
#2 <i>T. equigenitalis</i> (streptomycin resistant)	+ (Ct = 19,5)	-	-	-
#3 <i>T. equigenitalis</i> (streptomycin sensitive) + <i>S. aureus</i>	+ (Ct = 19,9)	-	-	-
#4 <i>S. aureus</i>	-	-	-	-
#5 <i>S. aureus</i>	-	-	-	-
#6 <i>K. pneumoniae</i> (Type 2) + <i>E. coli</i>	-	-	+ (Ct = 21)	-
#7 <i>T. equigenitalis</i> (streptomycin sensitive)	+ (Ct = 18,4)	-	-	-
#8 <i>T. equigenitalis</i> (streptomycin resistant) + <i>E. coli</i>	+ (Ct = 19,9)	-	-	-
#9 <i>E. coli</i>	-	-	-	-

PERFECT
CORRELATION
BETWEEN CULTURE
AND PCR RESULTS



Bio-X Diagnostics is ISO 9001:2008
certified to assure the best to its customers