

## **ADIAPURE™ Lysis FLEX**

**DIRECT LYSIS EXTRACTION KIT  
FOR THE DETECTION OF NUCLEIC ACIDS  
BY REAL-TIME ENZYMATIC GENE AMPLIFICATION (PCR TEST)**

**USE WITH ADIALYO RANGE ONLY**

**References:**  
ADPLF1-500



# ADIAPURE™ Lysis FLEX

|                                                              |   |
|--------------------------------------------------------------|---|
| MAIN CHANGE SINCE PREVIOUS VERSION .....                     | 3 |
| I. GENERAL INFORMATION .....                                 | 4 |
| 1. Purpose of the kit .....                                  | 4 |
| 2. Description of test .....                                 | 4 |
| II. MATERIAL & REAGENTS .....                                | 5 |
| 1. Composition of kit .....                                  | 5 |
| 2. Validity and storage .....                                | 5 |
| 3. Equipment required, but not supplied in the kit .....     | 5 |
| III. USE OF THE SAMPLES AND THE CONTROLS .....               | 6 |
| 1. Precautions .....                                         | 6 |
| 2. Storage of nucleic acid extracts .....                    | 6 |
| 3. Controls preparation .....                                | 6 |
| A. <i>Negative control of extraction (required)</i> .....    | 6 |
| B. <i>Positive control of extraction (recommended)</i> ..... | 6 |
| IV. ANALYSIS PROTOCOL .....                                  | 7 |
| 1. Extraction from swab .....                                | 7 |
| 2. Extraction from FTA CARD .....                            | 7 |
| 3. From faces and environmental sample (dung scraping) ..... | 8 |
| 4. Amplification .....                                       | 8 |
| V. INDEX OF SYMBOLS .....                                    | 9 |

## Main change since previous version

---

N/A                                      Not Applicable (first publication)  
Correction                                Correction of document anomalies  
Technical change                        Addition, revision and/or removal of information related to the product  
Administrative                            Implementation of non-technical changes noticeable to the user  
Note: minor typographical, grammar and formatting changes are not included in the revision historic

| Release Date | Part Number | Change type | Change summary                                                                              |
|--------------|-------------|-------------|---------------------------------------------------------------------------------------------|
| 2020/05      | NELF1-01    | N/A         | First publication                                                                           |
| 2023/06      | NELF1-02    | Technical   | Addition of Influenza virus and Infectious Bronchitis Virus detection from swabs            |
| 2023/10      | NELF1-03    | Technical   | Addition protocol from FTA CARD                                                             |
| 2024/03      | NELF1-04    | Technical   | Addition of <i>Mycobacterium avium</i> subsp. <i>paratuberculosis</i> detection from faeces |

## I. General information

---

### 1. Purpose of the kit

ADIAPURE™ Lysis FLEX is a nucleic acids extraction kit based on chemical lysis. The nucleic acids can be used without purification with the PCR amplification kit of **ADIALYO™** range validated with the kit.

### 2. Description of test

Adiogene validated the ADIAPURE™ Lysis FLEX kit from swabs, FTA cards, faeces, and environmental samples for the detection of pathogens in combination with the **ADIALYO™** range. The following table summarises the validated protocols.

|                           |                                                           | Swabs | FTA CARDS | Faeces | Environmental samples (dung scraping) |
|---------------------------|-----------------------------------------------------------|-------|-----------|--------|---------------------------------------|
| <b>Avian diseases</b>     | <i>Mycoplasma gallisepticum</i>                           | X     | X         |        |                                       |
|                           | <i>Mycoplasma synoviae</i>                                | X     | X         |        |                                       |
|                           | <i>Mycoplasma meleagridis</i>                             | X     | X         |        |                                       |
|                           | <i>Mycoplasma iowae</i>                                   | X     | X         |        |                                       |
|                           | Influenza Virus                                           | X     | X         |        |                                       |
|                           | Infectious Bronchitis Virus                               | X     | X         |        |                                       |
| <b>Ruminants diseases</b> | <i>Mycobacterium avium</i> subsp. <i>Paratuberculosis</i> |       |           | X      | X                                     |

For the pool size, refer to the user manual of the ADIALYO™ kit concerned.

## II. Material & reagents

---

### 1. Composition of kit

The ADIAPURE Lysis FLEX kit contains the following buffers:

| REF | ADIADP01S1-500 |              |                           |
|-----|----------------|--------------|---------------------------|
| LF1 | .....          | Lysis buffer | 5 x 100 mL (ready-to-use) |
| L2  | .....          | Enzyme       | 1 x 1.1 mL (ready-to-use) |
| LF3 | .....          | Lysis buffer | 1 x 28 mL (ready-to-use)  |

---

### 2. Validity and storage

On receipt, the kit should be aliquoted and stored at +2/8 °C. For better stability, it is recommended to keep L2 reagent at a temperature below -15 °C.

In this case, the kit is stable at least 1 year.

Do not mix reagents of two different batches.

### 3. Equipment required, but not supplied in the kit

**Warning: The material should be Nuclease-free (e.g. autoclaved 25 minutes twice at +120 °C or once 60 minutes at +121 °C).**

- Class II Microbiological Safety Cabinet.
- Incubator, heating bath or block heater.
- Vortex.
- 1 - 10 µL pipette, 20 - 200 µL pipette and 200 - 1000 µL pipette.
- Nuclease-free filter tips.
- Nuclease-free microtubes: 1.5 mL and 2 mL.
- Sterile tubes of 5, 10 or 15 mL.
- Powder-free latex or nitrile gloves.
- Fecal preparation system:
  - ADIAPREP™ (Bio-X Diagnostics ref.: 200 tests, ADPREP-200).
- Grinding beads.
  - For Mixer Mill:
    - ADIAPURE™ ALIQUOTED GLASS BEADS (Bio-X Diagnostics, 480 tests: ref. ADIADPBIA-480).
    - ADIAPURE GLASS BEADS RACKS 4x96 (Bio-X Diagnostics, 384 tests: ref ADPBIA-4x96).
  - For Rycoblyser and Fast prep:
    - Lysing Matrix B (MP biomedical, 100 tubes, ref. 116911100).

### III. Use of the samples and the controls

---

#### 1. Precautions

**Caution:**

Prepare the buffers of the kit according to the §II.2.

The buffers could contain toxic substances, please consult the MSDS safety data sheet.

The storage temperature must be respected.

We strongly recommend that only appropriately trained personnel perform this extraction. Ensure the accuracy and precision of the micropipettes used. The quality of the obtained results depends upon rigorous respect of good laboratory practices.

The PCR generates large amount of amplified DNA. A few molecules of amplified products are sufficient to generate a positive result. Do not open the PCR tubes after amplification.

Samples for analysis should be handled and disposed of as biological waste. **Take all measures of security and confinement required for the manipulation of the concerned biological agents.**

Before starting the process, read the entire protocol and scrupulously respect it.

#### 2. Storage of nucleic acid extracts

Extracted nucleic acids are temperature sensitive molecules. Crude extracts should be stored at the end of extraction on melting ice or at +2/8 °C for up to 24 hours, then at <-15 °C.

#### 3. Controls preparation

Several controls should be included per trial of analysis.

The use of different controls included in the ADIALYO™ kits allows validating all the analysis process steps (extraction and amplification) for all the samples.

- The endogenous or exogenous internal control included in the ADIALYO™ kits allow validation of the extraction and amplification steps for every sample.
- The positive control included in the ADIALYO™ kits allows validation of the amplification of the specific target.

Other controls should or must be added:

##### A. Negative control of extraction (required)

To verify the absence of cross-contamination, at least one negative control must be included per trial (e.g. AFNOR NF U47-600-1 guidelines suggest including a negative control per 24 columns centrifuged or four negative controls per trial of 96-wells plate). The control is a negative sample, for example a buffer used for dilution.

##### B. Positive control of extraction (recommended)

A positive control could be added in each trial. The control is a sample including the specific pathogen. It could come from a positive sample available in the laboratory or from a negative sample spiked with a solution of the specific pathogen. This positive control will be closed to the limit of detection of the method. It will inform about the fidelity of the obtained results between different trials.

## IV. Analysis protocol



Mix the buffers before use.

### 1. Extraction from swab

| Samples     | Mycoplasmas / Influenza Virus / Infectious Bronchitis Virus                                                                                                                                                                                                                                                                                                                                                                                                      |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Preparation | <p>Cut <b>1 to 6 swabs</b> in a 5 mL tube</p> <p>Add <b>1 mL of LF1 buffer</b> in the case of 1 to 3 swabs analysis<br/><b>Or add 2 mL of LF1 buffer</b> in the case of 4 to 6 swabs analysis</p> <p>Mix by vortexing 10 sec / tube</p> <p>Transfer <b>50 µL</b> of the supernatant in a microtube or in a PCR microplate well</p>                                                                                                                               |
| Lysis       | <p>Add <b>50 µL of LF3 buffer</b> and <b>2 µL of L2 buffer*</b></p> <p>Cover and homogenize</p> <p>Incubate <b>5 minutes at +65 °C</b> then <b>15 minutes at +95 °C</b></p> <p>Allow cooling to ensure the accuracy of subsequent pipetting.</p> <p><i>Note: Add <b>0,5 µL of EPC-Ext/PCR well</b> during the amplification step for the concerned kits.<br/>Refer to the corresponding package insert for preparation, storage, and use of the control.</i></p> |

\*Prior to use, a pre-mix of both reagents can be prepared and dispensed into each sample.

EPC-Ext is provided in the kit for the pathogen of interest.

Refer to the corresponding package insert for the preparation, storage, and use of the control.

### 2. Extraction from FTA CARD

| Samples | Mycoplasmas / Influenza Virus / Infectious Bronchitis Virus                                                                                                                                                                                                                                                                                                                               |
|---------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Lysis   | <p>Cut a 3 mm<sup>2</sup> piece of the FTA Card</p> <p>Place it in a microtube</p> <p>Add 100 µL of LF1 buffer, vortex 10 sec and discard the liquid</p> <p>Add 100 µL of LF3 buffer, vortex 10 sec and discard the liquid</p> <p>Add 100 µL of LF3 buffer, vortex 10 sec and incubate <b>5 minutes at +95 °C</b></p> <p>Allow cooling to ensure the accuracy of subsequent pipetting</p> |

### 3. Extraction from faeces and environmental sample (dung scraping)

| Sample      | <i>Mycobacterium avium</i> subsp. <i>paratuberculosis</i>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Preparation | <p>Collect <b>1 spoonful</b> of faecal matter using the ADIAPREP device. Scrape off the spoonful and reinsert it into the device. Vortex until a homogeneous suspension is obtained.</p> <p>Transfer <b>1 mL</b> into a microtube or 96-well plate.</p> <p>Centrifuge for <b>5 minutes at 3000 g</b> and remove the supernatant with a pipette.</p> <p>Add 300 mg of grinding beads and <b>500 µL of LF1 buffer</b> to the pellet.</p> <p>Grind for <b>5 minutes at 30 Hz</b> on Mixer Mill or 3 x 45 seconds at 4 m/sec on Fast Prep/Ribolyser then <b>centrifuge for 5 minutes at 3000 g</b>.</p> <p>Collect <b>50 µL</b> of supernatant.</p> |
| Lysis       | <p>Add <b>50 µL of LF3 buffer, 2 µL of L2 buffer and 5 µL of EPC-Ext*</b></p> <p>Cover and homogenize</p> <p>Incubate <b>5 minutes at +65 °C</b> then <b>15 minutes at +95 °C</b></p> <p>Allow cooling to ensure the accuracy of subsequent pipetting</p>                                                                                                                                                                                                                                                                                                                                                                                       |

\*Prior to use, a pre-mix of both reagents can be prepared and dispensed into each sample.

EPC-Ext is provided in the kit for the pathogen of interest.

Refer to the corresponding package insert for the preparation, storage, and use of the control.

### 4. Amplification

For the amplification of extracted nucleic acids, please refer to “Amplification Protocol” and “Reading and Interpretation” paragraphs of the **ADIALYO™** user manual of the pathogen of interest.

## V. Index of symbols

---

| Symbol                                                                            | Meaning                           |
|-----------------------------------------------------------------------------------|-----------------------------------|
|  | Catalog number                    |
|  | Manufacturer                      |
|  | Upper temperature limit           |
|  | Use by date                       |
|  | Batch code                        |
|  | Consult Instructions for Use      |
|  | Contains sufficient for <n> tests |
|  | Keep away from sunlight           |
|  | For veterinary in vitro use only  |

Bio-X Diagnostics, the logos, ADIAGENE, ADIALYO™, ADIAPURE™ and ADIAVET™ are used, pending and/or registered trademarks belonging to ADIAGENE and/or Bio-X Diagnostics, or one of its subsidiaries, or one of its companies. Any other name or trademark is the property of its respective owner.



**ADIAGENE S.A.S.**  
9, rue Gabriel Calloët-Kerbrat  
22440 Ploufragan - France

RCS 417 876 299  
Tel. +33 (0)2 96 68 40 20  
[www.biox.com](http://www.biox.com)