



ADIAVET™ TOXO FAST TIME

TEST FOR THE DETECTION OF *TOXOPLASMA GONDII* BY REAL-TIME
ENZYMATIC GENE AMPLIFICATION (PCR TEST)

References:

ADI273-50 (50 reactions)
ADI273-100 (100 reactions)



ADIAVET™ TOXO FAST TIME

I.	REVISION HISTORIC.....	3
II.	GENERAL INFORMATIONS	4
1.	Purpose of the test.....	4
2.	<i>Toxoplasma gondii</i>	4
3.	Description and purpose of the test.....	4
III.	MATERIAL & REAGENTS.....	5
1.	Reagents provided with the kit.....	5
2.	Validity and storage.....	5
3.	Use of controls	5
	A. <i>Use of EPC-Ext</i>	5
	B. <i>Use of TOXO CTL +</i>	5
4.	Equipment required but not supplied in the kit.....	5
IV.	RECOMMENDATIONS BEFORE THE ANALYSIS OF SAMPLES.....	7
1.	Precautions.....	7
2.	Storage of samples and DNA extracts.....	7
3.	Samples preparation	7
	A. <i>From swab</i>	7
	B. <i>From tissue</i>	7
	C. <i>From brain</i>	8
	D. <i>From fecal sample</i>	8
4.	Controls preparation	8
V.	EXTRACTION AND PURIFICATION	9
1.	Extraction using QIAamp® DNA Mini kit	9
2.	Extraction using NucleoSpin® Tissue kit	10
3.	Extraction using ADIAMAG - DNA/RNA magnetic beads kit	11
VI.	AMPLIFICATION	12
VII.	INTERPRETATION OF RESULTS	13
1.	Definitions.....	13
2.	Validation and interpretation of results.....	13
	A. <i>Validation of the run</i>	13
	B. <i>Result interpretation</i>	14
VIII.	INDEX OF SYMBOLS	15

I. Revision historic

N/A	Not Applicable (first publication)
Correction	Correction of document anomalies
Technical change	Addition, revision and/or removal of information related to the product
Administrative	Implementation of non-technical changes noticeable to the user

Note: minor typographical, grammar and formatting changes are not included in the revision historic

Release Date	Part Number	Change type	Change summaru
2019-05	NE273-01	N/A	First publication

II. General informations

1. Purpose of the test

ADIAVET™ TOXO FAST TIME kit is intended to detect *Toxoplasma gondii* parasite using real-time Polymerase Chain Reaction (PCR) technology from swab, tissue (brain, placenta, heart, liver, spleen, kidney, lungs of fetus) specimens of bovine, ovine and caprine and from cat faeces.

2. *Toxoplasma gondii*

Toxoplasma gondii is a widespread protozoan parasite that can determine serious disease in human, mammals and birds. *Toxoplasma gondii* may cause foetal resorption, abortion at any stage of pregnancy, foetal mummification and stillbirth. In non-pregnant animals, majority of infections are asymptomatic or of little consequence.

T. gondii is a cyst-forming coccidian parasite. Two asexual stages of development occur in the intermediate host (cattle, horse), tachyzoite and bradyzoite, and one sexual stage (oocyst) occurs in the definitive host (cat).

Routine diagnosis of *T. gondii* abortion is made by observation of clinical signs (abortion), by histology of foetal tissue and by indirect fluorescent antibody test (IFAT), and ELISA tests on maternal and foetal blood. Today, PCR based methods are useful diagnostic tools commonly used for detection of the parasite in aborted foetus.

3. Description and purpose of the test

This test is based on enzymatic gene amplification or PCR technology.

Amplified products are detected in real-time thanks to a specific labelled hydrolysis probe (5'-exonuclease technology).

The ADIAVET™ TOXO FAST TIME kit enables the simultaneous detection of:

- *Toxoplasma gondii* (probe labelled in FAM),
- GAPDH, an internal control of extraction and amplification steps specific from an endogenous DNA (probe labelled with a fluorochrome read in the same spectra as VIC and HEX).

or

- EPC-Ext: an exogenous control for extraction and amplification added to cat feces samples.

Adiagene recommends using this test with DNA purification kits (Bio-X Diagnostics, Qiagen and Macherey-Nagel). Other purification kits can be used if they have been validated by the user.

Analysis options according to the specimen:

Specimen	Individual analysis
Swab (vaginal, placental...)	☒
Tissue (placenta, foetal tissues...)	☒
Brain	☒
Cat faeces	☒

III. Material & reagents

1. Reagents provided with the kit

REF ADI273-50		
A5	Amplification Solution	1 x 500 µl tube green caps (ready to use)
TOXO CTL+	Positive control <i>Toxoplasma gondii</i>	1 tube purple caps (dehydrated)
EPC-Ext	external control of extraction for fecal samples	1 x 300 µl tube yellow caps (ready to use)
NF-Water	Nuclease-free Water	1 x 1000 µl tube with blank caps (Ready to use)
REF ADI273-100		
A5	Amplification Solution	2 x 500 µl tubes green caps (ready to use)
TOXO CTL+	Positive control <i>Toxoplasma gondii</i>	1 tube purple caps (dehydrated)
EPC-Ext	external control of extraction for fecal samples	2 x 300 µl tubes yellow caps (ready to use)
NF-Water	Nuclease-free Water	1 x 1000 µl tube with blank caps (Ready to use)

Package insert downloadable from www.biox.com

2. Validity and storage

On receipt, the kit should be stored at <-15°C.

It is recommended to make fractions of A5 solution if it should be defrosted more than 3 times.

Do not defrost reagents more than 3 times.

Realtime reagents are susceptible to light: store them in the darkness.

The A5 reagent is ready to use for PCR reaction.

Do not mix reagents of two different batches.

3. Use of controls

A. Use of EPC-Ext

EPC-Ext is an external control of extraction used to analysis from fecal sample.

Make fractions of EPC-Ext and store them at <-15°C. It should not be defrosted more than 3 times.

For each faeces, it is recommended to add **5 µl** of **EPC-Ext** per sample.

B. Use of TOXO CTL+

TOXO CTL+ is a positive control of amplification.

Add **200 µl** of **NF-water** to the **TOXO CTL+** tube and vortex at least 20 seconds.

Aliquot this solution by 6 or 12 µl and store them to <-15°C.

For each analysis, use **5 µl** of **TOXO CTL+** in one of the wells.

4. Equipment required but not supplied in the kit

Material should be Nuclease-free (e.g. autoclaved 25 minutes twice at +120°C or once 60 minutes at +121°C)

- Thermal cycler with consumables for real-time PCR: 0.2 ml PCR tubes or closed 96-wells PCR plates with optical quality.
- Class II Microbiological Safety Cabinet
- Centrifuge for microtubes
- Etuve, heating baths or block heaters
- Vortex
- Grinder (Mixer Mill or Fast Prep)
- 1 - 10 µl pipette, 20 - 200 µl pipette and 200 - 1000 µl pipette
- Nuclease-free filter tips
- Nuclease-free microtubes: 1.5 ml and 2 ml

- Sterile tube of 5, 10 or 15 ml
- Latex or nitrile powder-free gloves
- Scalpel blades
- 96-100% ethanol solution
- Nuclease-free water
- PBS buffer
- Sterile saline water (NaCl 8.5 g/L)
- MEM solution
- Grinding beads:
 - ADIAPURE™ ALIQUOTED GLASS BEADS (Bio-X Diagnostics, ref. ADIADPBIA-192 (192 tests), ref. ADIADPBIA-480 (480 tests)), only for disruption equipment such as Mixer Mill
 - Lysing Matrix B tubes (MP Biomedical, 100 extractions: ref. 116911.100) only for grinder Fast Prep.
- **DNA extraction kit (individual silica columns)**
 - QIAamp® DNA Mini Kit (Qiagen, 50 tests: ref. 51304 or 250 tests: ref. 51306)
 - NucleoSpin® Tissue (Macherey-Nagel, 50 tests: ref.740952.50 or 250 tests: ref. 740952.250)

or

- **Automated DNA/RNA extraction kit using magnetic beads**
 - ADIAMAG (Bio-X Diagnostics, 200 tests: ref.NADI003)

IV. Recommendations before the analysis of samples

Before starting the test, read the entire protocol and scrupulously respect it.

1. Precautions

Adiogene has elaborated this PCR test with the use of Bio-X Diagnostics, Qiagen and Macherey-Nagel extraction kits. Other extraction kits can be used with a previous validation.

Follow the supplier's instructions for the storage, the preparation and the use of the extraction reagents.

Some kits include and/or need the use of toxic reagents. These reagents should be handled with gloves and into a chemical cabinet.

We strongly recommend that only appropriately trained personnel perform this test. Ensure the accuracy and precision of the micropipettes used. The quality of the obtained results depends upon rigorous respect of good laboratory practices.

The PCR generates large amount of amplified DNA. A few molecules of amplified products are sufficient to generate a positive result. It is important to reserve 2 rooms, one for manipulation of samples to be tested, and another one for amplified products analysis. Do not open the PCR tubes after amplification.

Samples for analysis should be handled and disposed of as biological waste. **Take all measures of security and confinement required for the manipulation of the concerned biological agents.**

We recommend using fractions of demineralised and saline water and to autoclave them 25 minutes at +120°C. Take a new fraction for each new manipulation to avoid contamination.

2. Storage of samples and DNA extracts

Samples can be stored a couple of days at +2/8°C. After 2 days, we recommend to store them at <-15°C.

Extracted DNAs are quite sensitive molecules. Crude extracts should be stored at the end of extraction on melting ice or at +2/8°C for 24 hours, then at <-15°C.

3. Samples preparation

A. From swab

Mix the **swab** with **1 ml** of **1X PBS buffer**.

Transfer **200 µl** in a previously identified microtube.

See § IV for the extraction and purification of DNA.

B. From tissue

Particular case of placentas:

May contain a large amount of microorganisms, manipulate them with extreme precaution.

1st method

Cut the cotyledon with a scalpel, then rub inside with a swab.

Perform analysis according to swab protocol.

2nd method

Perform analysis according to tissue protocol.

Place **20-30 mg** of minced tissues in a previously identified microtube.
See § IV for the extraction and purification of DNA.

C. From brain

Put **one volume** of brain and **one volume** of MEM or sterile saline water. Homogenise.
When more than 10g of brain sample is analysed, it is possible to put the sample in a bag for blender (Mix 2 or Blender with paddles)
Transfer **200 µl** supernatant in a previously identified microtube.
See § IV for the extraction and purification of DNA.

D. From fecal sample

Place **1 g** of faeces in a previously identified 10 ml or 15 ml sterile tube.
Add **10 ml** of PBS 1X (*This preparation is stable 24h at room temperature*).
Vortex until a homogenous solution is obtained
Allow to settle 2 to 5 minutes.
Transfer **500 µl** of supernatant in a previously identified microtube.
Centrifuge at 3 000 g for 5 minutes. Discard the supernatant.
Homogenise the pellet with **1 ml** of PBS 1X (*This solution is stable 24 hours at room temperature*).
Transfer **500 µl** of the supernatant in a containing 300 mg of glass beads
Disrupt 10 minutes at 30 Hz with a Mixer Mill (or transfer the obtained solution in a Lysis Matrix B tube and disrupt 3 x 45 seconds at 4 m/sec with the Fast Prep).
Centrifuge 5 minutes at 15 000 g.
Transfer **200 µl** supernatant in a previously identified microtube.
Add **5 µl** of **EPC-Ext** in the sample.
See § IV for the extraction and purification of DNA.

4. Controls preparation

The use of controls allows verifying the reliability of the results.
The controls are included per trial of analysis. A trial is defined as all the samples treated in the same conditions.

All the steps of the analysis procedure (extraction+amplification), for all the types of samples, are validated with the association of the controls included in the kit.

- The internal endogenous control (GAPDH) naturally found in the samples verifies the extraction and amplification steps of each sample (except fecal sample).
- The external control (EPC-Ext) verifies the extraction and amplification steps of fecal sample.
- The TOXO CTL+ validates the amplification of the target.

Other controls must or could be added:

- **Negative control of extraction (required)**

To verify the absence of cross-contamination, at least one negative control must be included per trial. The control is a negative sample, for example a buffer used for dilutions.

- **Positive control of extraction (recommended)**

A positive control could be added in each trial. The control is a sample including *T. gondii*. It could come from a positive sample available in the laboratory or from a negative sample spiked with a solution of *T. gondii*. This positive control will be closed to the limit of detection of the method. It will inform about the fidelity of the obtained results between different trials.

V. Extraction and Purification

1. Extraction using QIAamp® DNA Mini kit

All the centrifugations are performed at room temperature.

Preparation of the sample	Take 200 µl of liquid sample or 20/30 mg of minced tissues , prepared as described before.
Lysis	Add 180 µl of ATL buffer , 20 µl of proteinase K . Vortex. Incubate 30 minutes at +70°C (or a night at +56°C).
	Add 200 µl of AL buffer . Vortex. Incubate 10 minutes at +70°C .
Binding preparation	Add 200 µl of ethanol 100% . Homogenise the mixture by pipeting (~10 times) or by vortex (~15 secondes).
Transfer to columns and binding to the membrane	Identify columns, apply the whole obtained solution to the corresponding column. Centrifuge 1 minute at 10 000 g. <i>If the whole sample has not been loaded once, apply the residual volume onto the column and centrifuge 1 minute at 10 000 g.</i>
1 st wash	Change the collection tube and add 500 µl of AW1 buffer to the column. Centrifuge 1 minute at 10 000 g.
2 nd wash	Change the collection tube and add 500 µl of AW2 buffer to the column. Centrifuge 1 minute at 10 000 g.
Column dry step	Change the collection tube. Centrifuge 3 minutes at 10 000 g.
Elution	Transfer the column to a microtube. Add 200 µl of AE buffer . Incubate ~1 minute at room temperature and centrifuge 1 minute at 10 000 g.
Storage	Close the tubes, identify and store at +2/8°C for 24 hours, then at <-15°C .

2. Extraction using NucleoSpin® Tissue kit

All the centrifugations are performed at room temperature.

Preparation of the sample	Take 200 µl of liquid sample or 20/30 mg of minced tissues, prepared as described before.
Lysis	Add 180 µl of T1 buffer , 25 µl of proteinase K . Vortex. Incubate 30 minutes at +70°C (or a night at +56°C).
	Add 200 µl of B3 buffer . Vortex. Incubate 10 minutes at +70°C .
Binding preparation	Add 200 µl of ethanol 100% . Homogenise the mixture by pipeting (~10 times) or by vortexing (~15 secondes).
Transfer to columns and binding to the membrane	Identify columns, apply the whole obtained solution to the corresponding column. Centrifuge at 10 000 g/1 minute . <i>If the whole sample has not been loaded once, apply the residual volume onto the column and centrifuge 1 minute at 10 000 g.</i>
1 st wash	Change the collection tube and add 500 µl of BW buffer to the column. Centrifuge 1 minute at 10 000 g .
2 nd wash	Change the collection tube and add 600 µl of B5 buffer to the column. Centrifuge 1 minute at 10 000 g .
Column dry step	Change the collection tube. Centrifuge 3 minutes at 10 000 g .
Elution	Transfer the column to a microtube. Add 200 µl of BE buffer . Incubate ~1 minute at room temperature and centrifuge 1 minute at 10 000 g .
Storage	Close the tubes, identify and store at +2/8°C for 24 hours , then at <-15°C .

3. Extraction using ADIAMAG - DNA/RNA magnetic beads kit

See the NEKF user manual available of the web site mentioned on the certificate of analysis included in the used ADIAVET™ kit.

VI. Amplification

a - Determine the number of analysed samples including the controls (e.g. positive and negative extraction controls, positive control of amplification (CTL+) and No Template Control (NTC)).

b - Defrost the A5 solution at room temperature. Vortex. Dispense **10 µl** of the A5 solution in each PCR tubes or PCR plate wells with a micropipette with a Nuclease-free tip.

c- **Immediately replace the A5 tube at <-15°C and in darkness.**

d- For each sample, the extraction negative control (obligatory) and the extraction positive control (recommended) add **5 µl** of purified extract to the **10 µl** of A5 solution.

For the CTL+, add **5 µl** of the solution obtained in § II-3 to the **10 µl** of A5 solution.

For the No Template Control (NTC), nothing is added to the A5 solution.

Immediately replace purified DNA extracts at +2/8°C or at <-15°C. Take care to have no bubbles in the bottom of the wells.

e- Once all the tubes have been prepared, run the real-time PCR amplification.

The *Toxoplasma gondii* target is read in FAM. The Internal Control is read in VIC or HEX. The Quencher is non fluorescent. The solution contains a passive reference ROX for the ABI machines. Fluorescence is read during the elongation step (1 minute at 60°C).

The following program is defined for **ABI Prism** thermalcyclers (like 7500, StepOne...) from **Applied Biosystems** (check the "emulation 9600" option if available), for the **MX3000P**, **MX3005P** and **AriaMx** from **Agilent**, for the **LightCycler 480** from **Roche**, for the **Rotorgene** from **Qiagen**, and for the **CFX96-Touch Real-Time PCR Detection System** from **Biorad**:

Standard program		Short program	
2 minutes at 50°C		2 minutes at 95°C	
10 minutes at 95°C			
15 seconds at 95°C**	45 cycles	5 seconds at 95°C	45 cycles
1 minute at 60°C		30 seconds at 60°C *	

* 32 seconds for the ABI7500 Applied Biosystems

** 30 seconds for the MX3005P

Roche diagnostic: LightCycler 2, LightCycler 480 :

NOTE: The use of *LightCycler* thermalcyclers requires a calibration manipulation. *Adiagene* will furnish process chart and reagents required for this calibration.

Contact us if you wish to use other thermalcyclers.

VII. Interpretation of results

1. Definitions



The « **base line** » corresponds to the background of fluorescence and qualifies the non characteristic part of the curve observed during the first cycles.

The « **Characteristic amplification curve** » qualifies a fluorescence curve with an exponential phase, a linear phase and a plateau phase.

The « **threshold line** » has to be placed over the background in the exponential phase of a characteristic amplification curve or a group of characteristic amplification curves.

The « **threshold cycle** » (**Ct**) of a well corresponds, for each detected fluorophore, at the crossing point of the threshold line with the fluorescent curve. The Ct value expressed by the machine for each well depends on the threshold position and on the quantity of target sequences initially present in the PCR well.

Example of characteristic amplification curve



2. Validation and interpretation of results

Display the FAM curves from the plate and set the threshold value as indicated above.

Proceed in the same mean for the HEX (or in the same spectra VIC) curves.

A. Validation of the run

Amplification is **valid** if the following results are obtained for the controls:

Controls	No Template Control (NTC)	Amplification positive control (CTL+)	Extraction negative control	Extraction positive control *
FAM amplification	No	yes	no	yes
HEX amplification	No	no/yes	no	no/yes
Validation of	Absence of contamination for amplification	Amplification of the target	Absence of contamination for extraction	Extraction and amplification steps

* Optional

The indicative Ct values (FAM and HEX dyes) of the positive control (CTL+) are indicated in the certificate of analysis of the kit.

B. Result interpretation

DNA extraction and amplification for each sample are **valid** if at least a characteristic amplification curve is observed for *Toxoplasma gondii* (FAM) or for the internal control (HEX).

Example	A	B	C	D
FAM amplification	no	yes	yes	no
HEX amplification	yes	no	yes	no
Result	negative	positive	positive	undetermined

The sample is considered as **negative** if a characteristic amplification curve is observed in HEX without any amplification in FAM (example A).

The sample is considered as **positive** if a characteristic amplification curve is observed in FAM (example B). Internal control can be co-amplified (example C).

A total absence of characteristic amplification curve for a sample (example D) shows a defective DNA extraction (lost or destruction of DNA) or a deficient real-time PCR (inhibitors in the sample, program error or no template added). In this case, we recommend first to repeat the test with pure and tenfold diluted DNA in Nuclease-free water. Then, if the test is still not valid, a new extraction is recommended.

VIII. Index of symbols

Symbol	Meaning
	Catalogue number
	Manufacturer
	Upper temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Keep away from sunlight
	For veterinary in vitro use only – For animal use only

Bio-X Diagnostics, the logos, ADIAGENE and ADIAVET™ are used, pending and/or registered trademarks belonging to ADIAGENE and/or Bio-X Diagnostics, or one of its subsidiaries, or one of its companies. Any other name or trademark is the property of its respective owner.



 **S.A.R.L. ADIAGENE**
9, rue Gabriel Calloët-Kerbrat
22440 Ploufragan - France

RCS 417 876 299
Tel. +33 (0)2 96 68 40 20
www.biox.com