

## ADIAPURE™ Lysis FLEX

**DIRECTE LYSIS EXTRACTION KIT  
FOR THE DETECTION OF NUCEIC ACIDS  
BY REAL-TIME ENZYMATIC GENE AMPLIFICATION (PCR TEST)**

**USING ONLY WITH ADILYO RANGE**

**References:**

ADPLF1-500



# ADIAPURE™ Lysis FLEX

<b>MAIN CHANGE SINCE PREVIOUS VERSION .....</b>	<b>3</b>
<b>I. GENERAL INFORMATION .....</b>	<b>4</b>
1. Purpose of the kit.....	4
2. Description of test.....	4
<b>II. MATERIAL &amp; REAGENTS .....</b>	<b>5</b>
1. Composition of kit .....	5
2. Validity and storage.....	5
3. Equipment required, but not supplied in the kit.....	5
<b>III. USE OF THE SAMPLES AND THE CONTROLS .....</b>	<b>6</b>
1. Precautions.....	6
2. Storage of nucleic acid extracts.....	6
3. Controls preparation .....	6
A. <i>Negative control of extraction (required)</i> .....	6
B. <i>Positive control of extraction (recommended)</i> .....	6
<b>IV. ANALYSIS PROTOCOL.....</b>	<b>7</b>
1. Extraction from swab.....	7
2. Amplification.....	7
<b>V. INDEX OF SYMBOLS.....</b>	<b>8</b>

## Main change since previous version

---

N/A                                      Not Applicable (first publication)  
Correction                                Correction of document anomalies  
Technical change                        Addition, revision and/or removal of information related to the product  
Administrative                            Implementation of non-technical changes noticeable to the user  
Note: minor typographical, grammar and formatting changes are not included in the revision historic

Release Date	Part Number	Change type	Change summary
2020/05	NELF1-01	N/A	First publication

## I. General information

---

### 1. Purpose of the kit

ADIAPURE™ Lysis FLEX is a nucleic acids DNA extraction kit based on the chemical lysis. The nucleic acids can be used without purification with the PCR amplification kit of **ADILYO™** range validated with the kit.

### 2. Description of test

Adiogene validated the ADIAPURE™ Lysis FLEX kit from swabs for the detection of DNA pathogens in combination with the **ADILYO™** range.

The following table summarises the validated protocols.

		Swab
<b>Avian diseases</b>	<i>Mycoplasma gallisepticum</i>	X
	<i>Mycoplasma synoviae</i>	X
	<i>Mycoplasma meleagridis</i>	X
	<i>Mycoplasma iowae</i>	X

For the pool size, refer to the user manual of the ADILYO™ kit concerned.

## II. Material & reagents

---

### 1. Composition of kit

The ADIAPURE Lysis FLEX kit contains the following buffers:

---

<b>REF</b> ADIADP01S1-500		
LF1 .....	Lysis buffer	5 x 100 mL (ready-to-use)
L2 .....	Enzyme	1 x 1.1 mL (ready-to-use)
LF3 .....	Lysis buffer	1 x 28 mL (ready-to-use)

---

### 2. Validity and storage

On receipt, the kit should be aliquoted and stored at +2/8°C. For better stability, it is recommended to keep reagent L2 at a temperature below -15°C.

In this case, the kit is stable at least 1 year.

Do not mix reagents of two different batches.

### 3. Equipment required, but not supplied in the kit

**Warning: The material should be Nuclease-free (e.g. autoclaved 25 minutes twice at +120°C or once 60 minutes at +121°C)**

- Class II Microbiological Safety Cabinet
- Incubator, heating bath or block heater
- Vortex
- 1 - 10 µL pipette, 20 - 200 µL pipette and 200 - 1000 µL pipette
- Nuclease-free filter tips
- Nuclease-free microtubes: 1.5 mL and 2 mL
- Sterile tubes of 5, 10 or 15 mL
- Powder-free latex or nitrile gloves

### III. Use of the samples and the controls

---

#### 1. Precautions

**Caution:**

Prepare the buffers of kit according to the §II.2.

The buffers could contain toxic substances, please consult the MSDS safety data sheet.

The storage temperature must be respected.

We strongly recommend that only appropriately trained personnel perform this extraction. Ensure the accuracy and precision of the micropipettes used. The quality of the obtained results depends upon rigorous respect of good laboratory practices.

The PCR generates large amount of amplified DNA. A few molecules of amplified products are sufficient to generate a positive result. Do not open the PCR tubes after amplification.

Samples for analysis should be handled and disposed of as biological waste. **Take all measures of security and confinement required for the manipulation of the concerned biological agents.**

Before starting the process, read the entire protocol and scrupulously respect it.

#### 2. Storage of nucleic acid extracts

Extracted DNAs are quite sensitive molecules. Crude extracts should be stored at the end of extraction on melting ice or at +2/8°C for max. 24 hours, then at <-15°C.

#### 3. Controls preparation

Several controls should be included per trial of analysis.

The mix of the different controls included in the Adiavet kits allows validating all the steps (extraction and amplification) of the analysis process for all the samples.

- The endogenous or exogenous internal control included in the ADIAVET™ kits allows validating the extraction and amplification steps of each sample.
- The positive control included in the ADIAVET™ kits allows validating the amplification of the specific target.

Other controls should or must be added:

##### A. Negative control of extraction (required)

To verify the absence of cross-contamination, at least one negative control must be included per trial (e.g. AFNOR NF U47-600-1 guidelines suggest to include a negative control per 24 columns centrifuged or four negative controls per trial of 96-wells plate). The control is a negative sample, for example a buffer used for dilution.

##### B. Positive control of extraction (recommended)

A positif control could be added in each trial. The control is a sample including the specific pathogen. It could come from a positive sample available in the laboratory or from a negative sample spiked with a solution of the specific pathogen. This positive control will be closed to the limit of detection of the method. It will inform about the fidelity of the obtained results between different trials.

## IV. Analysis protocol

---

### 1. Extraction from swab

 Mix the buffers before use.

Samples	Mycoplasmas
<b>Preparation</b>	<p>Cut <b>1 to 6 swabs</b> in a 5 mL tube.</p> <p>Add <b>1 mL of LF1 buffer</b> if 1 to 3 swabs analyse <b>Or add 2 mL of LF1 buffer</b> if 4 to 6 swabs analyse</p> <p>Mix by vortexing 10 sec / tube</p> <p>Transfer <b>50 µL</b> of the supernatant in a microtube or in a well of PCR microplate</p>
<b>Lysis</b>	<p>Add <b>50 µL of LF3 buffer</b> and <b>2 µL of L2 buffer*</b></p> <p>Cover and homogenize.</p> <p>Incube <b>5 minutes at +65°C</b> then <b>15 minutes at +95°C</b>.</p> <p>Let to cool, to ensure the accuracy of subsequent pipetting.</p>

\*Just before using, a pre-mix of the both reagents could be prepared then added to each sample.

### 2. Amplification

For the amplification of extracted nucleic acids, please refer to “Amplification Protocol” and “Reading and Interpretation” paragraphs of the user manual **ADILYO™** of the pathogen of interest.

## V. Index of symbols

---

Symbol	Meaning
	Catalogue number
	Manufacturer
	Upper temperature limit
	Use by date
	Batch code
	Consult Instructions for Use
	Contains sufficient for <n> tests
	Keep away from sunlight
	For veterinary in vitro use only

Bio-X Diagnostics, the logos, ADIAGENE, ADIAPURE™ and ADIAVET™ are used, pending and/or registered trademarks belonging to ADIAGENE and/or Bio-X Diagnostics, or one of its subsidiaries, or one of its companies. Any other name or trademark is the property of its respective owner.

---



**ADIAGENE S.A.S.**  
9, rue Gabriel Calloët-Kerbrat  
22440 Ploufragan - France

RCS 417 876 299  
Tel. +33 (0)2 96 68 40 20  
[www.biox.com](http://www.biox.com)