

MONOSCREEN[®] AbELISA

User manual

 BIO K 451_NO_(EN)_V02
31/08/2022

MonoScreen AbELISA *Neospora caninum* EASY

Reference : BIO K 451

ELISA Kit for the serodiagnosis of bovine neosporosis

Monowell, indirect test

For *in vitro* and strictly veterinary use

Sample	Species
Blood serum	Bovine
Individual milk (skimmed* and non-skimmed)	Bovine

* 20 min. 4000 g centrifugation

To order

Product reference	BIO K 451/5
Format	5 plates, 8-well strip
Reactions	480 tests

Kit composition

	BIO K 451/5
Microplate	5
Washing solution (20X)	1 X 250 ml
Colored dilution buffer (1X)	2 X 100 ml
Conjugate (50X)	1 X 1,4 ml
Reference - Positive serum - serum (black cap)	1 X 0,5 ml
Reference - Positive serum - milk (yellow cap)	1 X 0,5 ml
Reference - Negative serum	1 X 0,5 ml
Reference - Tracer	1 X 0,5 ml
Single component TMB solution (1X)	1 X 55 ml
Stopping solution (1X)	1 X 30 ml

Revision history

n/a	first version
V1.1.0	addition of the «individual milk» matrix
V1.2.0	amendment to item «G. Validation of results»
V02	removal of item «BIO K 451/2»

Note : minor changes to typography, grammar and formatting are not included in the revision history.


 Bio-X Diagnostics is ISO 9001
certified to assure the best to its customers

 Bio-X Diagnostics - 38, Rue de la Caestienne (PAE) - 5580 Rochefort - Belgium
 Tel : 0032(0)84.32.23.77 - Fax : 0032(0)84.31.52.63
 www.biox.com - E-mail : info@biox.com


A. Introduction

Neospora caninum is a protozoan initially described as a parasite of the dog in which it is responsible for myositis and encephalitis. Bovine neosporosis is now recognized as a major cause of abortion in cattle. It is strongly suspected in 20% of farms with repeated abortion and a seropositive cow for *Neospora caninum* is 3 times more likely to have an abortion than a seronegative cow. Vertical transmission is standard (at least 80% of calves from seropositive cows are contaminated).

B. Test principle

96-well microplates were sensitized by a specific monoclonal antibody of a *Neospora caninum* protein. The antibody ensures the capture and purification of this protein from a protozoan lysate.

Blood sera and milks are diluted in the dilution buffer. After incubation and washing of the preparation, the conjugate is added, a specific monoclonal antibody coupled with peroxidase. At the end of a second incubation of 30 minutes at $21\pm 3^{\circ}\text{C}$ and a second wash, the revelation solution is added (single component TMB solution). If specific immunoglobulins anti-*Neospora caninum* are present in the serum or milk, the conjugate remains attached to the well containing the protozoan and the enzyme catalyzes the transformation of colorless chromogene into a blue product. The intensity of the colouring is proportional to the specific antibody content in the sample.



C. Additional material and required equipment (not provided)

- Distilled/demineralized water
- Graduated mono- or multichannel pipettes (2-20µl, 20-200µl and 100-1000µl range) and single-use tips
- Microplate reader (450nm filter)
- Microplate washer and shaker (optional)
- Dilution microplate
- Standard laboratory equipment : graduated cylinder, tube rack, lid, ...

D. Precautions for use

- The reagents must be kept between $+2$ and $+8^{\circ}\text{C}$. The wash solution may be stored at room temperature.
- Unused strips must be stored with the desiccant in the hermetically sealed aluminium envelope.
- Do not use reagents beyond shelf-life date.
- Do not use reagents from other kits.
- Make sure to use distilled/demineralized water.
- The stopping solution contains 1 M phosphoric acid. Handle it carefully.
- Used material must be disposed of in compliance with the legislation in force regarding environmental protection and biological waste management.
- Keep the TMB solution away from light.

E. Preparation of solutions

- The washing solution must be diluted 20-fold in distilled/demineralized water. The cold solution crystallizes spontaneously. Bring the vial to $21\pm 3^{\circ}\text{C}$ to make sure that all crystals have disappeared; mix the solution well and withdraw the necessary volume.
- The dilution buffer is ready to use. The dilution buffer is coloured in yellow.
- The conjugate must be diluted 50-fold in the dilution buffer.
- The stopping solution is ready to use.
- La TMB solution is ready to use. It must be perfectly colourless.

F. Procedure

- Bring all the reagents to $21\pm 3^{\circ}\text{C}$ before use.
- Carefully read through the previous points.

N.B. : To avoid differences in incubation time between samples, sample dilutions and reference dilutions can be prepared in a dilution microplate before transfer (200 µl) into the test microplate using a multi-channel pipette.

Serum protocol (1/20 dilution)

- 1. Distribute 190 µl/well of dilution buffer. Add 10 µl of serum and references per well. Homogenize by pipetting up and down.
- 2. Cover and incubate the plate at 21 ± 3°C during 30 ± 3 min.
- 3. Remove the content of the microplate. **Wash the microplate 3 times with 300 µl of washing solution per well.** Avoid the formation of bubbles in the wells and the dessication of the microplate between each wash.
- 4. Add **100 µl of diluted conjugate** per well. Cover with a lid and incubate the plate at 21 ± 3°C during 30 ± 3 min.

Milk protocol (1/4 dilution)

- 1. For milk (1/4 dilution): distribute the dilution solution at rate of 150 µl per well. Add samples at rate of 50 µl per well. Homogenize by pipetting up and down.

For the references (1/20 dilution): distribute 190 µl of dilution buffer per well. Add 10µl per well of reference. Homogenize by pipetting up and down.
- 2. Cover and incubate the plate at 21 ± 3°C during 60 ± 5 min.
- 3. Remove the content of the microplate. **Wash the microplate 3 times with 300 µl of washing solution per well.** Avoid the formation of bubbles in the wells and the dessication of the microplate between each wash.
- 4. Add **100 µl of diluted conjugate** per well. Cover with a lid and incubate the plate at 21 ± 3°C during 60 ± 5 min.

Joint protocol

- 5. Remove the content of the microplate. **Wash the microplate 3 times with 300 µl of washing solution per well.** Avoid the formation of bubbles in the wells and the dessication of the microplate between each wash.
- 6. Distribute **100 µl of TMB solution** per well.
- 7. Incubate at 21 ± 3°C during 10 ± 1 min away from the light, without covering.
- 8. Distribute the stopping solution at rate of **50 µl per well.** Colour changes from blue to yellow.
- 9. Record the optical densities using a plate spectrophotometer with a 450 nm filter within 5 minutes after adding the stopping solution.

G. Validation of results

The test can only be validated if :

- the difference between positive and negative serum optical density readings is greather than 0,450.
$$OD_{\text{positive serum (serum or milk)}} - OD_{\text{negative serum}} > 0,450$$
- the negative serum gives an optical density of less than 0,400.
$$OD_{\text{negative serum}} < 0,400$$

H. Interpretation of results

Calculate for each sample its coefficient (S/P %) using the following formula:

$$S/P \% = \frac{OD_{\text{sample}} - OD_{\text{negative serum}}}{OD_{\text{positive serum (serum or milk)}} - OD_{\text{negative serum}}} * 100$$

Results		Status
Serum	S/P % < 70 %	Negative
	S/P % ≥ 70 %	Positive
Milk	S/P % < 50 %	Negative
	S/P % ≥ 50 %	Positive

Get the interpretation of your results quickly and easily using **AnalysisScreen**, our free online platform, available on our website : <https://www.biox.com>

ANALYSISCREEN


- AnalysisScreen™ is the new module for reading and interpreting all types of **Monoscreen™** and **Multiscreen™** ELISA plates. AnalysisScreen™ is :
- Free
 - Accessible online via our website : <https://www.biox.com>
 - Updated in real time
 - Compatible with all Bio-X Diagnostics plate designs
 - Very easy to use



SCAN ME


Serum protocol

190 μ l of dilution buffer + 10 μ l of serum (1/20)
 190 μ l of dilution buffer + 10 μ l of references (1/20)



Milk protocol

150 μ l of dilution buffer + 50 μ l of milk (1/4)
 190 μ l of dilution buffer + 10 μ l of references (1/20)



Add 100 μ l of diluted conjugate (1/50 dilution)



Add 100 μ l of diluted conjugate (1/50 dilution)



Add 100 μ l of TMB solution



Add 50 μ l of stopping solution

Record the optical densities

