



BTv

Reference: ADL35Y1-100 & ADL35Y1-500

Test for the detection of Bluetongue virus (BTV) by real time enzymatic amplification

PCR Test – 100 & 500 reactions

For veterinary *in vitro* use only



Sample	Individual analysis	Pool of sample possible *, up to:
Blood	✓	20
Tissue	✓	✗

* Depending on the epidemiological case and on the quality of samples

Kit composition

Content		ADL35Y1-100 Kit	ADL35Y1-500 Kit
		100 reactions	500 reactions
A6	Amplification solution	1 lyophilized vial with white caps (To reconstitute)	5 lyophilized vials with white caps (To reconstitute)
Rehydration buffer	A6 rehydration solution	1 x 6 mL vial (Ready to use)	1 x 6 mL vial (Ready to use)
BTV CTL+	BTV positive control	1 tube with purple cap (To reconstitute)	2 tubes with purple cap (To reconstitute)
PCR Buffer	CTL+ rehydration solution	1 x 1000 µL tube with white cap (Ready to use)	1 x 1000 µL tube with white cap (Ready to use)

Revision history

Date	Version	Modifications
01/2023	V01	Adaptation of packaging in 100 reactions
02/2025	V02	Modification of "NF water" by "PCR Buffer". Addition of cervids and camelids species. Adaptation of packaging in 500 reactions.

Note: minor typographical, grammar and formatting changes are not included in the revision history.

A. Introduction

The bluetongue virus is a non-contagious viral arthropod-borne infectious disease due to an Orbivirus (family Reoviridae, virus ARN), mainly transmitted by hematophagous midges from *Culicoides* genus.

The clinical expression is widely dependent on the environmental parameters (nutritional state, parasitism and bacterial infections concomitant) and on the individual sensitivity. 27 distinct serotypes exist inducing partial or no cross protections between them. Transmission by pregnant ewes has also been described. Transmission by contaminated blood injection is possible when needles and syringes are re-used.

Samples for virus detection are bloods of animals with anticoagulants (EDTA). Virus is detected by isolation on embryonated eggs, *in vitro* cell culture, immunofluorescence on cell culture or by PCR.

B. Test principle

ADIALYO™ BTV test is based on the reverse transcription (RT) of RNA into complementary DNA. Then, cDNA is amplified with a DNA polymerase using specific primers of Bluetongue Virus. Both enzymatic reactions occur in the same tube (One-step RT-PCR). This test is intended to detect simultaneously, in one well:

- Bluetongue Virus (FAM labelled probe).
- GAPDH internal control of extraction and amplification specific from an endogenous nucleic acid (HEX labelled probe or its equivalent).

C. Storage conditions

- Store the kit at a temperature below +2/8 °C after reception.
- Store away from sunlight and keep dry.
- After reconstitution, prepare aliquots and store them at a temperature below -15 °C until the expiration date.
- Do not thaw more than 3 times.

D. Material required but not provided

- Real-time Thermal cycler and device.
- Instrument for homogenous mixing of tubes.
- Pipettes of 1 - 10 µL, 20 - 200 µL and 200 - 1000 µL.
- Nuclease-free filtered pipette tips.
- Nuclease-free microtubes of 1,5 mL and 2 mL.
- Powdered-free latex or nitrile gloves.
- Nuclease-free water.
- Kit for nucleic acids extraction.

Additional kits for method adoption and PCR

- **Extraction Positive Control BTV (Ref.: ADC35EPC).**
Supplier reference material for method adoption that can also be used as a sentinel (Calibrated between 1 and 100xLOD_{Method}).
- **LD_{PCR} Positive Control – BTV (Ref.: ADC35YLD)**
Confirmation of performances – LOD_{PCR} of kit.

E. Warnings and precautions

- For veterinary *in vitro* use only.
- For animal use only.
- For professional use only.
- All instructions must be read before performing the test and strictly respected.
- Do not use reagents after the expiration date.
- Do not use reagents if the packaging is damaged.
- Do not open PCR wells or tubes after amplification.
- Do not mix reagents from different batches.

- Used material must be disposed of in compliance with the legislation in force regarding environmental protection and biological waste management.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).

F. Nucleic acids extraction

1. Extraction kits

Nucleic acids must be extracted from the samples before using the kit. The RNA/DNA purification kits listed below are recommended by Bio-X Diagnostics:

Product name	Extraction system	Number of tests and reference
ADIAMAG	Magnetic beads	200 tests: ref. NADI003 800 tests: ref. NADI003-XL

For the extraction, consult the user manual version, available on the website, indicated on the certificate of analysis included in the PCR kit of interest.

Extraction protocols are described in validation data. Other purification kits can be used if they have been validated by the user.

After extraction, nucleic acid extracts can be kept on ice or at +2/8 °C for until use. For long term storage, they must be kept at a temperature below -15 °C or -65 °C.

2. Controls

Using controls allow to verify the reliability of the results. Controls can be included.

Control	Validation of	Usage
No Template Control (NTC)	Absence of amplification contamination	5 µL NF-Water in a well per run
BTB CTL+	BTB target amplification	5 µL CTL+ in a well per run
Negative extraction control	Absence of contamination for the extraction and amplification	1 extraction (water or lysis buffer) per run
Positive extraction control	Extraction and amplification	1 extraction (Positive sample between 1 et 100X LOD _{Method}) per run

G. Procedure

1. Amplification solution A6 preparation

- Add **1000 µL** of « **Rehydration buffer** » per A6 tube.
- Homogenize tube contents using a mixer, such as vortex, at least 20 seconds.
- After reconstitution, aliquot the solution and store at a temperature below -15 °C until the expiration date. Do not thaw more than 3 times.
- To use the A6, please refer to §« Amplification », Step 1.

2. Preparation of CTL+

- Add **200 µL** of « **PCR Buffer** » per tube.
- Homogenize the tubes by suction and pressure, then use a mixer, such as vortex, > 20 seconds and until dissolution of the blue pellet.
- After reconstitution, aliquot and store the solution at a temperature below -15 °C until the kit expiration date. Do not thaw more than 3 times.

For each assay, use **5 µL** of denatured CTL+ (see § « Denaturation of nucleic acid ») in one of the dedicated wells (see § « Amplification », Step 2).

3. Denaturation of nucleic acids extracts

- For each sample and control(s), transfer 10 µL of nucleic acids extracts into a tube or 96-plate and store the rest at a temperature below -15 °C or -65 °C.
- Incubate 3 minutes at +95 °C in a thermal cycler or heating block.
- Immediately transfer the tubes or 96-plate on melting ice or refrigerated block until use (to prevent RNA renaturation).
- To use, please refer to § « Amplification », Step 2.

4. Amplification

Warning:

- Before starting, rehydrate or thaw reagents at room temperature in the dark.
- Homogenize all reagents and samples before use.
- Store reagents at a temperature below -15 °C after use.

Step 1: Dispense **10 µL** of amplification solution (A6) per well.

Step 2: Dispense **5 µL** of denatured nucleic acids extracts and **5 µL** of controls in each dedicated well.

Use PCR Buffer for the No Template Control (NTC).

Step 3: Cover the wells with an appropriate optical film or caps.

Step 4: Set up the amplification program.

The following program is defined for ABI Prism thermocyclers (like 7500, QuantStudio5, Step-one...) from Applied Biosystems, for Mx3000, Mx3005P and AriaMx from Agilent, for LightCycler from Roche Diagnostics, for Rotor-Gene Q from Qiagen, for CFX96 and Chromo 4 from Biorad and for MIC from BioMolecular System.

DNA/RNA Program	
10 min. 45 °C	
2 min. 95 °C	
5 sec. 95 °C	40 cycles
30 sec. 60 °C*	

*Reading and parameters for fluorescence acquisition:

Fluorochrome	Absorbance (nm)	Emission (nm)
FAM	494	520
HEX or equivalent	538	554
ROX	575	602

Note: The Quencher is non-fluorescent. The A6 solution contains a passive reference read in the same spectra as ROX for ABI machines.

For other thermal cycler instruments, please contact your sales representative or the customer relations department.

H. Reading and interpretation

Display all curves and position the threshold line for each fluorochrome.

1. Test validation

Amplification is valid if the following results are obtained.

Expected Ct (Threshold Cycle) values for the CTL+ are indicated on the certificate of analysis of the kit.

Controls	Amplification		Validation of
	FAM	HEX or equivalent	
No Template Control (NTC)	No	No	Absence of amplification contamination
BTV CTL+	Yes	Yes/No	Target amplification
Extraction negative control	No	No	Absence of extraction contamination
Extraction positive control	Yes	Yes/No	Extraction and amplification steps

2. Results interpretation

Nucleic acids extraction and amplification are **valid** for each sample if at least one typical amplification curve is observed in FAM and/or HEX or equivalent.

Amplification		Interpretation
FAM	HEX or equivalent	Bluetongue virus
No	Yes	Undetected
Yes	Yes	Detected
Yes	No	Detected
No	No	Undetermined

« **Undetermined** »: no characteristic amplification curve.

Possible causes:

Defective PCR due to inhibitors, set up error, absence of samples, degraded samples and/or issue with nucleic acids extraction (loss or destruction of nucleic acids).

Recommendations:

Set up a new PCR assay using pure nucleic acids extracts and 10x dilutions in Nuclease-free water;

If the assay is inconclusive, perform a new nucleic acids extraction.

Symbols

Symbole	Signification
	Catalog number
	Manufacturer
	Temperature limitation
	Use by
	Batch code
	Consult Instructions for Use
	Contain sufficient for "n" tests
	For veterinary <i>in vitro</i> use only – For animal use only
	Keep away from sunlight
	Keep dry

- 1 Extract nucleic acids with

Adia^X
Mag



Scan me to
discover Adiamag™


- 2 Add **1000 µL** of **Rehydration buffer** to the **A6** amplification solution




- 3 Distribute **10 µL** of **A6** amplification solution

- 4 Denature **10 µL of nucleic acids**, and each control **3 min at 95°C**.
Transfer immediately to ice or ice packs.


- 5 Distribute **5 µL** of denatured **nucleic acids**, denatured **CTL+** and **PCR buffer**.



- 6 Seal the wells



- 7 Start PCR analysis



*The notes do not replace the instructions for use of which they are a summary.