



## LAWSONIA

Reference: ADL02Y1-100

Test for the detection of *Lawsonia intracellularis* by real time enzymatic amplification

PCR Test – 100 reactions

For veterinary *in vitro* use only



| Sample                           | Individual analysis | Pool of sample possible*, up to: |
|----------------------------------|---------------------|----------------------------------|
| Swab (rectal...)                 | ✓                   | 3                                |
| Tissue (intestinal mucosa, etc.) | ✓                   | ✗                                |
| Faeces                           | ✓                   | 3                                |
| Environmental sample             | ✓                   | ✗                                |
| Oral fluid                       | ✓                   | ✗                                |

\* Depending on the epidemiological case and on the quality of samples

## Kit composition

| Content            |                                                  | ADL02Y1-100 Kit                                         |
|--------------------|--------------------------------------------------|---------------------------------------------------------|
|                    |                                                  | 100 reactions                                           |
| A6                 | Amplification solution                           | 1 lyophilized vial with blank caps<br>(To reconstitute) |
| Rehydration buffer | Rehydration solution                             | 1 x 6 mL vial<br>(Ready to use)                         |
| LAW CTL+           | <i>Lawsonia intracellularis</i> positive control | 1 tube with purple cap<br>(To reconstitute)             |
| NF-Water           | Nuclease-Free Water                              | 1 x 1000 µL tube with white cap<br>(Ready to use)       |

## Revision history

| Date    | Version | Modifications |
|---------|---------|---------------|
| 01/2025 | V01     | First version |

Note: minor typographical, grammar and formatting changes are not included in the revision history.

## A. Introduction

*Lawsonia intracellularis* is the ethological agent of Proliferative Enteropathy (PE), causing a disease commonly known as 'Ileitis' and having a major economic impact on the pig industry. *L. intracellularis* has two clinical presentations in pigs:

Porcine hemorrhagic enteropathy (PHE) is the acute, haemorrhagic and severe form, affecting animals aged between 4 and 12 months. Porcine intestinal adenomatosis (PIA) is the most common form. It is a chronic infection diagnosed in pigs aged between 6 and 20 weeks. Infection by *L. intracellularis* is spread mainly via the faecal-oral route. Infected animals excrete the bacteria in their feces, and contamination often occurs in confined conditions, such as in intensive livestock farms.

## B. Test principle

ADIALYO™ LAWSONIA test is based on the amplification of specific *Lawsonia intracellularis* DNA. This test is intended to detect simultaneously, in one well:

- *L. intracellularis* (FAM labelled probe).
- Internal control of amplification specific from an exogenous DNA (HEX labelled probe or its equivalent).

## C. Storage conditions

- Store the kit at a temperature below +2/8 °C after reception.
- Store away from sunlight and keep dry.
- After reconstitution, prepare aliquots and store them at a temperature below -15 °C until the expiration date.
- Do not thaw more than 3 times.

## D. Material required but not provided

- Real-time Thermal cycler and device.
- Instrument for homogenous mixing of tubes.
- Pipettes of 1 - 10 µL, 20 - 200 µL and 200 - 1000 µL.
- Nuclease-free filtered pipette tips.
- Nuclease-free microtubes of 1,5 mL and 2 mL.
- Powdered-free latex or nitrile gloves.
- Nuclease-free water.
- Kit for nucleic acids extraction.

### Additional kits (on request) for method adoption and PCR

- **LD<sub>PCR</sub> Positive Control – Law (Ref.: ADC02LD)**  
Confirmation of performances – LOD<sub>PCR</sub> of kit.

## E. Warnings and precautions

- For veterinary *in vitro* use only.
- For animal use only.
- For professional use only.
- All instructions must be read before performing the test and strictly respected.
- Do not use reagents after the expiration date.
- Do not use reagents if the packaging is damaged.
- Do not open PCR wells or tubes after amplification.
- Do not mix reagents from different batches.
- Used material must be disposed of in compliance with the legislation in force regarding environmental protection and biological waste management.
- This kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not totally guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled observing the usual safety precautions (do not ingest or inhale).

## F. Nucleic acids extraction

### 1. Extraction kits

Nucleic acids must be extracted from the samples before using the kit. The RNA/DNA purification kits listed below are recommended by Bio-X Diagnostics:

| Product name        | Extraction system | Number of tests and reference                         |
|---------------------|-------------------|-------------------------------------------------------|
| ADIAMAG             | Magnetic beads    | 200 tests: ref. NADI003<br>800 tests: ref. NADI003-XL |
| ADIAPURE Lysis Flex | Direct lysis      | 500 mL: ref. ADPLF1-500                               |

For the extraction, consult the user manual version, available on the website, indicated on the certificate of analysis included in the PCR kit of interest.

Extraction protocols are described in validation data. Other purification kits can be used if they have been validated by the user.

After extraction, nucleic acid extracts can be kept on ice or at +2/8 °C for a few hours or until use. For long term storage, they must be kept at a temperature below -15 °C or -65 °C.

### 2. Controls

Using controls allows to verify the reliability of the results.

| Control                     | Validation of                                                 | Usage                                                                           |
|-----------------------------|---------------------------------------------------------------|---------------------------------------------------------------------------------|
| No Template Control (NTC)   | Absence of amplification contamination                        | 5 µL NF-Water in a well per run                                                 |
| LAW CTL+                    | LAW target amplification                                      | 5 µL CTL+ in a well per run                                                     |
| Negative extraction control | Absence of contamination for the extraction and amplification | 1 extraction (water or lysis buffer) per run                                    |
| Positive extraction control | Extraction and amplification                                  | 1 extraction (Positive sample between 1 et 100X LOD <sub>METHOD</sub> ) per run |

## G. Procedure

### 1. Amplification solution A6 preparation

- Add **1000 µL** of « **Rehydration buffer** » per A6 tube.
- Homogenize tube contents using a mixer, such as vortex, at least 20 seconds.
- After reconstitution, aliquot the solution and store at a temperature below -15 °C until the expiration date. Do not thaw more than 3 times.
- To use the A6, please refer to § « Amplification », Step 1.

### 2. Preparation of CTL+ control

- Add **200 µL** of « **NF-Water** » per tube.
- Homogenize the tubes by suction and pressure, then use a mixer, such as vortex, > 20 seconds and until dissolution of the blue pellet.
- After reconstitution, aliquot and store the solution at a temperature below -15 °C until the kit expiration date. Do not thaw more than 3 times.
- For each assay, use **5 µL** of CTL+ in one of the dedicated wells (see § « Amplification », Step 2).

### 3. Amplification

#### Warning:

- Before starting, rehydrate or thaw reagents at room temperature in the dark.
- Homogenize all reagents and samples before use.
- Store reagents at a temperature below -15 °C after use.

**Step 1:** Dispense **10 µL** of amplification solution (A6) per well.

**Step 2:** Dispense **5 µL** of nucleic acids extracts and **5 µL** of controls in each dedicated well.

Use NF-Water for the No Template Control (NTC).

**Step 3:** Cover the wells with an appropriate optical film or caps. Centrifuge briefly (recommended).

**Step 4:** Set up the amplification program.

The following program is defined for ABI Prism thermocyclers (like 7500, QuantStudio5, Step-one...) from Applied Biosystems, for Mx3000, Mx3005P and AriaMx from Agilent, for LightCycler from Roche Diagnostics, for Rotor-Gene Q from Qiagen, for CFX96 and Chromo 4 from Biorad and for MIC from BioMolecular System.

| DNA/RNA Program |           |
|-----------------|-----------|
| 10 min. 45 °C   |           |
| 2 min. 95 °C    |           |
| 5 sec. 95 °C    | 40 cycles |
| 30 sec. 60 °C*  |           |

\*Reading and parameters for fluorescence acquisition:

| Fluorochrome      | Absorbance (nm) | Emission (nm) |
|-------------------|-----------------|---------------|
| FAM               | 494             | 520           |
| HEX or equivalent | 530             | 549           |
| ROX               | 575             | 602           |

**Note:** The Quencher is non-fluorescent. The A6 solution contains a passive reference read in the same spectra as ROX for ABI machines.

For other thermal cycler instruments, please contact your sales representative or the customer relations department.

### H. Reading and interpretation

Display all curves and position the threshold line for each fluorochrome.

#### 1. Test validation

Amplification is valid if the following results are obtained.

Expected Ct (Threshold Cycle) values for the CTL+ are indicated on the certificate of analysis of the kit.

| Controls                    | Amplification |                   | Validation of                          |
|-----------------------------|---------------|-------------------|----------------------------------------|
|                             | FAM           | HEX or equivalent |                                        |
| No Template Control (NTC)   | No            | Yes               | Absence of amplification contamination |
| LAW CTL+                    | Yes           | Yes               | Target amplification                   |
| Extraction negative control | No            | Yes               | Absence of extraction contamination    |
| Extraction positive control | Yes           | Yes               | Extraction and amplification steps     |

### 2. Results interpretation

Nucleic acids extraction and amplification are **valid** for each sample if at least one typical amplification curve is observed in FAM and/or HEX or equivalent.

| Amplification |                   | Interpretation                  |
|---------------|-------------------|---------------------------------|
| FAM           | HEX or equivalent | <i>Lawsonia intracellularis</i> |
| No            | Yes               | Undetected                      |
| Yes           | Yes               | Detected                        |
| Yes           | No                | Detected                        |
| No            | No                | Undetermined                    |

« **Undetermined** » : no characteristic amplification curve.

#### Possible causes:

Defective PCR due to inhibitors, set up error, absence of samples, degraded samples and/or issue with nucleic acids extraction (loss or destruction of nucleic acids).

#### Recommendations:

Set up a new PCR assay using pure nucleic acids extracts and 10x dilutions in Nuclease-free water.

If the assay is inconclusive, perform a new nucleic acids extraction.

### Symbols

| Symbole | Signification                                                 |
|---------|---------------------------------------------------------------|
|         | Catalog number                                                |
|         | Manufacturer                                                  |
|         | Temperature limitation                                        |
|         | Use by                                                        |
|         | Batch code                                                    |
|         | Consult Instructions for Use                                  |
|         | Contain sufficient for "n" tests                              |
|         | For veterinary <i>in vitro</i> use only – For animal use only |
|         | Keep away from sunlight                                       |
|         | Keep dry                                                      |

1 Extract nucleic acids with

**Adia<sup>X</sup>**  
Pure



Scan me to  
discover Adiapure™

or

**Adia<sup>X</sup>**  
Mag



Scan me to  
discover Adiamag™


2 Add **1000 µL** of **Rehydration buffer** to  
the **A6** amplification solution






3 Distribute **10 µL** of **A6** amplification  
solution

4 Distribute **5 µL** of **nucleic acids**, **CTL+**  
and **NF-Water**

5 Seal the wells



6 Start PCR analysis



\*The notes do not replace the instructions  
for use of which they are a summary.